

FAMEs and CLA isomers



Column: **MEGA-10 - 0.25mm, 0.20 μ m, 100m**

Catalog Code: C-10-025-020-50-50


Conditions:

Oven Program: 40°C (4min), 10°C/min, 120°C (1min),
5°C/min, 180°C (18min), 2°C/min, 200°C (1min),
2°C/min, 240°C.

Carrier Gas: Hydrogen, 1.2mL/min, constant flow.

Injector: Split 250°C, 1:30 Split Ratio.


Detector: FID 250°C.



Analysis carried out with DANI Master GC.

MEGA[®]

www.mega.mi.it



Column: **MEGA-10 - 0.25mm, 0.20 μ m, 100m**

Catalog Code: C-10-025-020-50-50

Conditions:

Oven Program: 40°C (4min), 10°C/min, 120°C (1min), 5°C/min, 180°C (18min), 2°C/min, 200°C (1min), 2°C/min, 240°C.

Carrier Gas: Hydrogen, 1.2mL/min, constant flow.

Injector: Split 250°C, 1:30 Split Ratio.

Detector: FID 250°C.



contact us: info@mega.mi.it

 follow us [@MEGAColumns](https://twitter.com/MEGAColumns)

www.mega.mi.it

since
1980 